EL TESTRAIGHT STAIRLIFT

MADE IN USA Veteran Founded



BRING HOME **COMFORT & ACCESS**



ELITESTRAIGHT



COMBINE QUALITY, COMFORT & STYLE


Gain high performance, modern styling and heavy lift capacity with Bruno's topof-the line Elite straight stairlift. Glide up and down your stairs with a gentle, smooth ride from start to finish. Make your life a little easier with a Bruno Elite, featuring Made-In-America quality and maximum comfort.



The Bruno Elite's Limited Lifetime Warranty gives you the peace of mind to enjoy worry-free access to your entire home long into the future.

3 TOP CONSIDERATIONS

- ✓ Made in the USA Quality
- ✓ 400 lb. (190kg) Lift Capacity
- ✓ Limited Lifetime Warranty



HIGH-PERFORMANCE. EXCEPTIONAL RIDE.

Easy. Armrest control operates the stairlift up and down. Two wireless remote controls are also included.

Safe. Seat swivels 90° at top of steps and locks in place for safe exit away from stairs. Seat belt and obstruction sensors further ensure safety.

Comfortable. Generously sized seat is adjustable in height. Footrest height and space between armrests are also adjustable.

Space-Saving. Vertical rail leaves maximum open space on the steps. Fold arms, seat and footrest when not in use to save even more space.

Sleek. Covered gear rack on vertical rail gives a slim, clean appearance and reduces exposure to dirt and grease.

Dependable. Battery-powered performance, even during a power outage. Stairlift continually charges anywhere it stops on the rail.

Learn More bruno.com

\star \star \star \star \star \checkmark \checkmark Verified Customer

"Our Bruno stairlift experience was first class. From the initial contact, to the prompt well-done installation to the smooth operation, it was a +++ experience."

Ralph L.

MADE IN USA Veteran Founded

AVAILABLE OPTIONS

Add a manual or power folding rail when you have a doorway at the bottom of the stairs or an installation where the rail may cause an obstruction.

Get an extra-wide seat and/or larger footrest to best match specific body types.

Eliminate the need to lean over with a power folding footrest. Fold the seat, and the footrest will automatically fold, too. Change the standard swivel seat to power, and the seat rotates nearly 90 degrees at the top landing with a push of the control.

Request a four-point harness for additional securement points on the chair.

Multiple seat color options also available.










Not all options are shown.



ELTESTRAIGHT





TECH SPECS

Lift Capacity: 400 lb (181 kg)

Standard Seat Size: 18 in wide x 15 in deep (457 mm x 381 mm)

Larger Seat Option: 22.5 in wide x 16.5 in deep (572 mm x 419 mm)

Standard Footrest Size: 10 in long x 9 in wide (254 mm x 229 mm)

Larger Footrest Option: 12 in long x 11 in wide (305 mm x 279 mm)

Vertical Rail: Installs within 6 in (152 mm) of the wall

WARRANTY

Lifetime Warranty on Major Components -- motor, gearbox, drivetrain and stationary rail system -- for the lifetime of the Bruno Stairlift as long as the original purchaser owns the stairlift.

2-Year Warranty on all components after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

ABOUT BRUNO

Veteran-Founded | Family-Owned | Made in USA

Bruno makes stairlifts and vertical platform lifts to help people better access homes, public buildings and businesses, and scooter lifts to assist people in easily transporting their mobility devices.



© Bruno Independent Living Aids, Inc. 2022 2022-07-11-elite-mkt-sl-2.indd • *MKT-SL-2*